

=====  
 Acq. Operator : SYSTEM Seq. Line : 2  
 Acq. Instrument : LCMS-1 Location : P1-E3  
 Injection Date : 9/16/2020 7:34:09 AM Inj : 2  
 Inj Volume : 1.000 µl

Method : C:\CHEM32\1\METHODS\Purity195\_315.M  
 Last changed : 5/8/2020 11:20:44 AM by SYSTEM  
 Method Info : Operator: Oraganic QC lab  
 Gradient Method Starting @ 100 % Water/0 % Acetonitrile to 0 % Water/100 %  
 Acetonitrile in 13 min, then 3 min at 100/0, post-run time 2 min. Flow rate:0.3 ml  
 /min, column Zorbax Eclipse Plus C-18, 2.1x50 mm, 1.8 micron particule size.  
 Agilent 1290 Infinity II LC system, Agilent 6120 Quadropole MS Detector. No UV  
 Reference.

Additional Info : Peak(s) manually integrated

Current Chromatogram(s)

