

File :Q:\GCMS-08 Minimal\DATA\20190416\10- .D  
Operator : JP/KM  
Acquired : 16 Apr 2019 20:23 using AcqMethod KERLA-2017-100uL draw speed.M  
Instrument : GCMS 8  
Sample Name: M-8141A-1M  
Misc Info : Method 8141 Additions in 97.5/2.5 Hexane-Acetone  
Vial Number: 5



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

