

File :Q:\GCMS-08 Minimal\DATA\20190607\02- D
Operator : JP/KM
Acquired : 07 Jun 2019 17:17 using AcqMethod KERLA-2017_S50.M
Instrument : GCMS 8
Sample Name: M-8270-15
Misc Info : SV by Capillary Column GC/MS in 90:10 MeCl2/MeOH
Vial Number: 2



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

