

=====

Acq. Operator	: SYSTEM	Seq. Line	: 20
Acq. Instrument	: LCMS-1	Location	: P1-C1
Injection Date	: 6/16/2021 1:48:33 PM	Inj	: 1
		Inj Volume	: 1.000 µl
Acq. Method	: C:\CHEM32\1\METHODS\Purity195_315.M		
Last changed	: 5/21/2021 4:03:06 PM by SYSTEM		
Analysis Method	: C:\CHEM32\1\METHODS\Purity195_315.M		
Last changed	: 6/16/2021 10:33:12 AM by SYSTEM (modified after loading)		
Method Info	: Operator: Organic QC lab		
	Gradient Method Starting @ 100 % Water/0 % Acetonitrile to 0 % Water/100 % Acetonitrile in 13 min, then 3 min at 100/0, post-run time 2 min. Flow rate :0.3 ml/min, column Zorbax Eclipse Plus C-18, 2.1x50 mm, 1.8 micron particle size. Agilent 1290 Infinity II LC system, Agilent 6120 Quadropole MS Detector. No UV Reference.		

Additional Info : Peak(s) manually integrated

Current Chromatogram(s)

