

Sample Name : P-772S

=====

Acq. Operator : SYSTEM Seq. Line : 22  
Sample Operator : SYSTEM  
Acq. Instrument : LCMS-2 Location : P2-D-08  
Injection Date : 6/18/2021 4:24:13 PM Inj : 1  
Inj Volume : 1.000 µl



Acq. Method : C:\Users\Public\Documents\ChemStation\1\Data\06182101\06182101 2021-06-18  
09-46-17\Purity\_195-315.M

Last changed : 4/29/2021 10:41:07 AM by SYSTEM

Analysis Method : C:\Users\Public\Documents\ChemStation\1\Methods\Purity\_195-315.M

Last changed : 4/29/2021 10:41:07 AM by SYSTEM

Method Info : OPERATOR: Organic QC lab

General gradient method starting @100% water for 1 min @0.3 ml/min, the 100%  
Acetonitrile in 13 min, 100% Acetonitrile until T=16 min. Column Zorbax  
Eclipse plus C-18, 2.1mmx50mm, 1.8 µ particle size. Agilent 1290 Infinity  
II LC/MS system, Agilent 6125 quadrupole MS detector.

Additional Info : Peak(s) manually integrated

### Current Chromatogram(s)

