



=====

Acq. Operator : SYSTEM	Seq. Line : 8
Acq. Instrument : LCMS-1	Location : P1-D2
Injection Date : 6/6/2018 1:05:00 PM	Inj : 1
	Inj Volume : 1.000 µl

Quinmerac

Acq. Method : D:\DATA\06061801\06061801 2018-06-06 10-34-26\Purity195\_315.M  
Last changed : 6/6/2018 10:34:27 AM by SYSTEM  
Analysis Method : C:\Chem32\1\Methods\Purity195\_315.M  
Last changed : 6/6/2017 9:45:25 AM by SYSTEM  
Method Info : Operator: GA  
Gradient Method Starting @ 100 % Water/0 % Acetonitrile to 0 % Water /100 % Acetonitrile in 13 min, then 3 min at 100/0, post-run time 2 min . Flow rate:0.3 ml/min, column Zorbax Eclipse Plus C-18, 2.1x50 mm, 1.8 micron particule size. Agilent 1290 Infinity II LC system, Agilent 6120 Quadropole MS Detector. No UV Reference.

Additional Info : Peak(s) manually integrated

Current Chromatogram(s)

