



=====

Acq. Operator	: SYSTEM	Seq. Line	: 3
Acq. Instrument	: LCMS-1	Location	: P2-F3
Injection Date	: 7/23/2018 11:10:56 AM	Inj	: 1
		Inj Volume	: 1.000 µl

Method : C:\Chem32\1\Data\07231801 2018-07-23 10-31-58\PURITY195_315.M (Sequence Method)

Last changed : 7/23/2018 10:31:58 AM by SYSTEM

Method Info : Operator: GA
Gradient Method Starting @ 100 % Water/0 % Acetonitrile to 0 % Water/100 % Acetonitrile in 13 min, then 3 min at 100/0, post-run time 2 min. Flow rate:0.3 ml/min, column Zorbax Eclipse Plus C-18, 2.1x50 mm, 1.8 micron particule size. Agilent 1290 Infinity II LC system, Agilent 6120 Quadropole MS Detector. No UV Reference.

Current Chromatogram(s)

