

File :Q:\GCMS-08\DATA\20190214\35-  
Operator : JP/KM  
Acquired : 14 Feb 19 08:06 pm using AcqMethod KERLA-2017-100uL draw speed.M  
Instrument : GCMS 8  
Sample Name: CLP-031-R-WL  
Misc Info : CLP Surrogate Spiking Solution in MeOH  
Vial Number: 27



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

