

```
=====
Acq. Operator   : SYSTEM                               Seq. Line :    3
Acq. Instrument : LCMS-1                             Location  : P1-D2
Injection Date  : 5/12/2021 3:42:02 PM                Inj       :    1
                                                    Inj Volume: 1.000 µl
```



Different Inj Volume from Sample Entry! Actual Inj Volume : 2.000 µl

```
Method          : C:\CHEM32\1\METHODS\Purity195_315.M
Last changed    : 5/8/2020 11:20:44 AM by SYSTEM
Method Info     : Operator: Oraganic QC lab
                  Gradient Method Starting @ 100 % Water/0 % Acetonitrile to 0 % Water/100 %
                  Acetonitrile in 13 min, then 3 min at 100/0, post-run time 2 min. Flow rate
                  :0.3 ml/min, column Zorbax Eclipse Plus C-18, 2.1x50 mm, 1.8 micron
                  particule size. Agilent 1290 Infinity II LC system, Agilent 6120 Quadropole
                  MS Detector. No UV Reference.
```

Additional Info : Peak(s) manually integrated

Current Chromatogram(s)

