

File :D:\MassHunter\GCMS\1\DATA\07312301\19- ALR.D
Operator : Organic QC department
Acquired : 31 Jul 2023 16:04 using AcqMethod KERLA-2017.M
Instrument : 5977B MSD
Sample Name: M-525.2-TS-A
Misc Info : Method 525.2 - Tuning Standard in Acetone
Vial Number: 90



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

