



=====

Acq. Operator : SYSTEM	Seq. Line :
Sample Operator : SYSTEM	
Acq. Instrument : LCMS-2	Location : P1-F-10
Injection Date : 11/5/2021 3:47:45 PM	Inj : 2
	Inj Volume : 1.000 µl

Acq. Method : C:\Users\Public\Documents\ChemStation\1\Data\11052101\11032101 2021-11-05  
15-26-51\Purity\_195-315.M

Last changed : 4/29/2021 10:41:07 AM by SYSTEM

Analysis Method : C:\Users\Public\Documents\ChemStation\1\Methods\Purity\_195-315.M

Last changed : 4/29/2021 10:41:07 AM by SYSTEM

Method Info : OPERATOR: Organic QC lab

General gradient method starting @100% water for 1 min @0.3 ml/min, the 100% Acetonitrile in 13 min, 100% Acetonitrile until T=16 min. Column Zorbax Eclipse plus C-18, 2.1mmx50mm, 1.8 µm particle size. Agilent 1290 Infinity II LC/MS system, Agilent 6125 quadrupole MS detector.

Sample Info : 1,3,5-Trinitrobenzene @ 100ug/mL in AcCN/MeOH 50/50

Additional Info : Peak(s) manually integrated

Current Chromatogram(s)

