

=====

Acq. Operator	: SYSTEM	Seq. Line	: 7
Sample Operator	: SYSTEM		
Acq. Instrument	: LCMS-2	Location	: P1-F-07
Injection Date	: 12/1/2021 2:21:13 PM	Inj	: 1
		Inj Volume	: 1.000 µl

Acq. Method : C:\Users\Public\Documents\ChemStation\1\Data\12012101\12012101 2021-12-01 12-26-36  
Purity\_195-315.M

Last changed : 4/29/2021 10:41:07 AM by SYSTEM

Analysis Method : C:\USERS\PUBLIC\DOCUMENTS\CHEMSTATION\1\METHODS\Purity\_195-315.M

Last changed : 12/1/2021 8:27:34 AM by SYSTEM  
(modified after loading)

Method Info : OPERATOR: Organic QC lab  
General gradient method starting @100% water for 1 min @0.3 ml/min, the 100% Acetonitrile in 13 min, 100% Acetonitrile until T=16 min. Column Zorbax Eclipse plus C-18, 2.1mmx50mm, 1.8 µm particle size. Agilent 1290 Infinity II LC/MS system, Agilent 6125 quadrupole MS detector.

Additional Info : Peak(s) manually integrated

Current Chromatogram(s)

