



=====
 Acq. Operator : SYSTEM Seq. Line : 2
 Sample Operator : SYSTEM
 Acq. Instrument : LCMS-2 Location : P1-A11
 Injection Date : 7/22/2022 9:29:34 AM Inj : 2
 Inj Volume : 2.000 µl

Acq. Method : C:\Users\Public\Documents\ChemStation\1\Data\07222201\07212201 2022-07-22
 09-09-42\PFSA.M

Last changed : 1/26/2021 10:30:23 AM by SYSTEM

Analysis Method : C:\USERS\PUBLIC\DOCUMENTS\CHEMSTATION\1\METHODS\Purity_195-315.M

Last changed : 3/23/2022 11:04:54 AM by SYSTEM

Method Info : OPERATOR: Organic QC lab

General gradient method starting @100% water for 1 min @0.3 ml/min, the 100%
 Acetonitrile in 13 min, 100% Acetonitrile until T=16 min. Column Zorbax
 Eclipse plus C-18, 2.1mmx50mm, 1.8 µm particle size. Agilent 1290 Infinity
 II LC/MS system, Agilent 6125 quadrupole MS detector.

Sample Info : Perfluorooctane sulfonamide in MeOH

Current Chromatogram(s)

