

File : B:\QC Data\GCMS-08 Minimal\DATA\20190927\20-.D
Operator : ZA / KM / JP
Acquired : 27 Sep 2019 15:55 using AcqMethod KERLA-2017.M
Instrument : GCMS 8
Sample Name: ALR-048S-ET-10X
Misc Info :
Vial Number: 135

Column: DB-5, 30m, 0.25 ID, 0.25um
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

