

File : Q:\GCMS-06 Minimal\DATA\120121\61- YX.D
Operator : Organic QC Lab
Acquired : 03 Dec 2021 02:19 using AcqMethod CHICK_2019_2S_2uL.M
Instrument : GCMS 6
Sample Name: ALR-168S
Misc Info : Triethanolamine @100 ug/ml in MeOH
Vial Number: 122



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

