

File :Q:\GCMS-04 Minimal\DATA\021420A\110D
Operator :
Acquired : 18 Feb 20 9:48 am using AcqMethod VOC-100SPLIT-FAST.M
Instrument : GCMS4
Sample Name: ALR-CFC-001S-2X
Misc Info : RR
Vial Number: 31



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

