

File :Q:\GCMS-08 Minimal\DATA\20230905\11\_JP\_.D  
Operator : Organic QC Lab  
Acquired : 05 Sep 2023 16:48 using AcqMethod Hy\_PAH\_24cm\_S200.M  
Instrument : GCMS 8  
Sample Name: AOCs-006N  
Misc Info :in chloroform  
Vial Number: 9



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 2:1, 2µl inj.; GC/MS; injector 270°C

