

File : B:\QC Data\GCMS-04 Minimal\DATA\101019A\03-D
Operator : KM/ZA
Acquired : 10 Oct 2019 13:29 using AcqMethod VOC-50SPLIT DECKER.M
Instrument : GCMS4
Sample Name : APP-9-007-W
Misc Info : RS
Vial Number : 21



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

