

File : L:\QC Data\GCMS-04 Minimal\DATA\121819B\53-.D
Operator :
Acquired : 19 Dec 2019 12:03 using AcqMethod VOC- 12 SPLIT DECKER.M
Instrument : GCMS4
Sample Name : APP-9-015
Misc Info :
Vial Number : 76



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

