

File :Q:\GCMS-06 Minimal\DATA\102021\06- KM.D  
Operator : Organic QC Lab  
Acquired : 20 Oct 2021 16:11 using AcqMethod CHICK\_2019\_VOC-SV\_S100\_20min.M  
Instrument : GCMS 6  
Sample Name: APP-9-027  
Misc Info : bis(2-Chloroethyl)ether in MeOH  
Vial Number: 138



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

