


File : B:\QC Data\GCMS-08 Minimal\DATA\20190909\15-.D  
Operator : JJP / KM /  
Acquired : 09 Sep 2019 18:22 using AcqMethod KERLA-2017.M  
Instrument : GCMS 8  
Sample Name: APP-9-179  
Misc Info : PHENOL IN DCM A  
Vial Number: 103

 **AccuStandard**<sup>®</sup>  
Leader in Analytical Reference Standards  
Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 80c, 17c/min to 340c, 8min  
GC Parameters: Cons. Split, 12psi constant flow  
Split 20:1, 1uL inj.; GC/MS; INJ 270c

