

File : L:\QC Data\GCMS-08 Minimal\DATA\20220602\10 SAK.D
Operator : Organic QC Lab
Acquired : 02 Jun 2022 12:22 using AcqMethod Hy_Pest_15cm_S100.M
Instrument : GCMS 8
Sample Name: AS-E0011
Misc Info :
Vial Number: 126



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

