

File :Q:\GCMS-06 Minimal\DATA\091619\19-.D
Operator : D.C
Acquired : 16 Sep 19 09:34 pm using AcqMethod CHICK_2019_S100.M
Instrument : GCMS 6
Sample Name: AS-E0023
Misc Info : 1,2-Dichlorobenzene in MeOH
Vial Number: 37



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

