

File : Q:\GCMS-04 Minimal\DATA\112723\11- GMA.D  
Operator :  
Acquired : 27 Nov 2023 17:54 using AcqMethod VOC-100 SPLIT #H2 12MIN80DEG.M  
Instrument : GCMS4  
Sample Name: AS-E0043  
Misc Info :  
Vial Number: 6



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm  
Oven Program: 30°C, Hold 5min, 7°/m-60°, Hold 0min, 11°/m-225, Hold 1min  
GC Parameters: 1µL inj, GC-MS, Injector temp 250°C, Carrier Gas: Hydrogen

