

File :D:\MassHunter\GCMS\1\data\20230413B\72- -=JP.D
Operator : Organic QC Lab
Acquired : 14 Apr 23 06:43 pm using AcqMethod Hy_Pest_24cm_S10.M
Instrument : GCMS 8
Sample Name: AS-E0194
Misc Info : 2-Fluorobiphenyl; in acet
Vial Number: 51



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

