

File : Q:\GCMS-06 Minimal\DATA\062521\11- KM.D
Operator : Organic QC Lab
Acquired : 25 Jun 2021 13:48 using AcqMethod CHICK_2017_13MIN_S100.M
Instrument : GCMS 6
Sample Name: AS-E0250
Misc Info : o-Cresol in MeOH
Vial Number: 45



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

