

File : B:\QC Data\GCMS-04 Minimal\DATA\100819\66-.D
Operator : RM / KM/ZA
Acquired : 9 Oct 2019 21:42 using AcqMethod VOC-200SPLIT DECKER.M
Instrument : GCMS4
Sample Name : AS-E0305
Misc Info :
Vial Number: 3



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

