


File : L:\QC Data\GCMS-08 Minimal\DATA\20200320\12-.D  
Operator : Organic QC Lab  
Acquired : 20 Mar 2020 15:40 using AcqMethod Hy\_PAH\_24cm\_S100.M  
Instrument : GCMS 8  
Sample Name: AS-E0322  
Misc Info :  
Vial Number: 143

 **AccuStandard®**  
Leader in Analytical Reference Standards  
Column: DB-5MS, 30m, 0.25 ID, 0.25 um  
Oven Program: 80c, 17c/min to 340c, 8min  
GC Parameters: Cons. Split, 12psi constant flow  
Split 20:1, 1uL inj.; GC/MS; INJ 270c

