

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-06 Minimal\DATA\060222\0  
... 4- KM.D  
Operator : Organic QC Lab  
Instrument : GCMS 6  
Acquired : 02 Jun 2022 13:49 using AcqMethod VOC20Min#.M  
Sample Name: AS-E0360  
Misc Info : Carbon tetrachloride in MeOH



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

