

File : Q:\GCMS-08 Minimal\DATA\20200624\30 -DC.D  
Operator : Organic QC Lab  
Acquired : 24 Jun 20 08:12 pm using AcqMethod Hy\_M502\_S100.M  
Instrument : GCMS 8  
Sample Name : AS-E0368  
Misc Info : 1,2,3-Trichloropropane in Methanol  
Vial Number: 46



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

