

File :Q:\GCMS-07 Minimal\DATA\10291901\38-.D
Operator : KD & DC
Acquired : 30 Oct 2019 01:40 using AcqMethod KERLA-2017_S100.M
Instrument : 5977B MSD
Sample Name: AS-E0485
Misc Info : N-Nitrosomorpholine in MeOH
Vial Number: 132



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

