

Chrom Perfect Chromatogram Report



C:\CPData6\DATA\062519(1) 0001.RAW

ASM-177

1	Tetrachloro-m-xylene
2	a-BHC
3	g-BHC
4	Heptachlor
5	Endosulfan I
6	Dieldrin
7	Endrin
8	p,p'-DDD
9	p,p'-DDT
10	Methoxychlor
11	Decachlorobiphenyl

Column: DB-5, 30m, 0.25 ID, 0.25um
Oven Program: 80c, 0min, 20c/min to 320c, final 4min.
GC Parameters: Splitless; 2ul inj; ECD 310c; Inlet 250c;

