

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\020223\8
 ... 3- RM.D
 Operator :
 Instrument : GCMS4
 Acquired : 4 Feb 2023 3:17 using AcqMethod VOC-100 SPLIT #H2 28MIN(HOT).M
 Sample Name: ASTM-P-0080
 Misc Info :



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

