

File :Q:\GCMS-06 Minimal\DATA\100721\03- -DC.D
Operator : Organic QC Lab
Acquired : 07 Oct 21 09:54 am using AcqMethod CHICK_2019_VOC-SV_S100.M
Instrument : GCMS 6
Sample Name: ASTM-P-0140-QC
Misc Info : Daily QC Std
Vial Number: 44



Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 35c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 100:1, 1uL inj.; GC/MS; INJ 270c

