

File : Q:\GCMS-08 Minimal\DATA\20230808\22- _JP~.D
Operator : Organic QC Lab
Acquired : 08 Aug 23 08:07 pm using AcqMethod Hy_Pest_24cm_S10.M
Instrument : GCMS 8
Sample Name: BP-024S
Misc Info : 2,4-Dibromophenol; in tol
Vial Number: 26



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

