


File :Q:\GCMS-08 Minimal\DATA\20220606\15_JP_.D
Operator : Organic QC Lab
Acquired : 06 Jun 2022 14:17 using AcqMethod Hy_Pest_24cm_S10.M
Instrument : GCMS 8
Sample Name: C-044S
Misc Info : 22'35'-Tetra-PCB; in iso
Vial Number: 54

 **AccuStandard®**
Leader in Analytical Reference Standards
Column: DB-5MS, 30m, 0.25 ID, 0.25 um
Oven Program: 80c, 17c/min to 340c, 8min
GC Parameters: Cons. Split, 12psi constant flow
Split 20:1, 1uL inj.; GC/MS; INJ 270c

