

File :D:\MassHunter\GCMS\1\data\20220525\23- (JP).D  
Operator : Organic QC Lab  
Acquired : 25 May 22 04:19 pm using AcqMethod Hy\_Pest\_24cm\_S1.M  
Instrument : GCMS 8  
Sample Name: C-110S  
Misc Info : 233'4'6-PentaPCB;  
Vial Number: 63



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

