

File : Q:\GCMS-08 Minimal\DATA\20210610\08- [ JP ].D  
Operator : Organic QC Lab  
Acquired : 10 Jun 21 09:00 pm using AcqMethod Hy\_Pest\_15cm\_S10.M  
Instrument : GCMS 8  
Sample Name: C-155S-TP  
Misc Info : 2,2',4,4',6,6'-Hexa-PCB; in isooctane  
Vial Number: 82



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

