

File :Q:\GCMS-08 Minimal\DATA\20220422\32- =JP=.D  
Operator : Organic QC Lab  
Acquired : 22 Apr 22 08:05 pm using AcqMethod Hy\_Pest\_24cm\_S1.M  
Instrument : GCMS 8  
Sample Name: C-194S-TP (222041486)  
Misc Info : 22'33'44'55'-Octa-PCB; in iso  
Vial Number: 59



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

