

File : Q:\GCMS-08 Minimal\DATA\20220729\2 -jopo.D  
Operator : Organic QC Lab  
Acquired : 29 Jul 2022 15:31 using AcqMethod Hy\_Pest\_24cm\_S1.M  
Instrument : GCMS 8  
Sample Name: C-260S-TP  
Misc Info : Aroclor 1260; in iso  
Vial Number: 5



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

