

File :Q:\GCMS-08 Minimal\DATA\20231016\54- (JP).D  
Operator : Organic QC Lab  
Acquired : 17 Oct 23 07:01 am using AcqMethod Hy\_PAH\_24cm\_S200.M  
Instrument : GCMS 8  
Sample Name: CCME-PAH  
Misc Info : in dcm  
Vial Number: 17



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 80c, 1m, 15c/min to 340c, final 2min  
GC Parameters: Split 10:1, 1ul inj.; GC/MS; INJ 270c  
(Actual analytical conditions may vary from these general parameters)

