

File : Q:\GCMS-04 Minimal\DATA\121119B\11-.D
 Operator :
 Acquired : 12 Dec 2019 19:29 using AcqMethod VOC-100SPLIT DECKER.M
 Instrument : GCMS4
 Sample Name: CLP-022-LC
 Misc Info :
 Vial Number: 96

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
 Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
 GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

