

File :Q:\GCMS-08 Minimal\DATA\20220608\20- (JP).D  
Operator : Organic QC Lab  
Acquired : 08 Jun 22 02:39 pm using AcqMethod Hy\_PAH\_24cm\_S100.M  
Instrument : GCMS 8  
Sample Name: CLP-HC-SV-MIX4  
Misc Info : 2mg/ml in dcm  
Vial Number: 52



Column: DB-5, 30m, 0.25 ID, 0.25um  
Oven Program: 100c, 1m, 30c/min to 320c, final 2min.  
GC Parameters: Pulsed split 40psi 0.5m, 10psi constant flow  
Split 20:1, 1ul inj.; GC/MS; INJ 250c.

