

File :Q:\GCMS-08 Minimal\DATA\20231031\20- (JP).D  
Operator : Organic QC Lab  
Acquired : 31 Oct 23 05:54 pm using AcqMethod Hy\_VOC-SV\_24cm\_S100.M  
Instrument : GCMS 8  
Sample Name: CLP-PI-2.5X  
Misc Info : in meoh  
Vial Number: 17



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm  
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)  
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow  
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

