

File :Q:\GCMS-04 Minimal\DATA\052820A\49--RM.D
Operator : ORGANIC QC DEPT.
Acquired : 29 May 2020 18:33 using AcqMethod VOC-200SPLIT HOT.M
Instrument : GCMS4
Sample Name: D-3606-17-CM
Misc Info :
Vial Number: 75



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

