

File : \\SVR-FS-01\Organic QC\QC Data\GCMS-04 Minimal\DATA\051122\6  
... 1RM.D  
Operator :  
Instrument : GCMS4  
Acquired : 12 May 2022 3:24 using AcqMethod VOC-200SPLIT HOT260.M  
Sample Name: D-3606-QC-IS2-25ML  
Misc Info :



Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

