

File :Q:\GCMS-04 Minimal\DATA\050123\01-RM.D
Operator :
Acquired : 1 May 2023 12:23 using AcqMethod VOC24-100 SPLIT #H2 20MIN175DEG.M
Instrument : GCMS4
Sample Name: D-5441
Misc Info :
Vial Number: 11



Column: DB-624UI, 20m x 0.18mm ID x 1.0µm
Oven Program: 30, Hold 5min, 7-60, 0 min Hold 11-225, 1min Hold
GC Parameters: 1µl inj.; GC-MS; Injector temp 250°C carrier gas Hydrogen

