

File :K:\QC Data\GCMS-04 Minimal\DATA\051721\05- RM.D  
Operator :  
Acquired : 17 May 2021 12:37 using AcqMethod VOC-200SPLIT HOT260.M  
Instrument : GCMS4(221051308)

Column: DB-624 UI, 30m x 0.25mm ID x 1.4µm  
Oven Program: 35°C (hold 5min), 11°C/min to 60°C, 22°C/min to 230°C (hold 4min)  
GC Parameters: Split 100:1, 1µl inj.; GC/MS; Injector temp 240°C

Misc Info : E0  
Vial Number: 38

