

File :Q:\GCMS-08 Minimal\DATA\20210922\20- _JP=.D
Operator : Organic QC Lab
Acquired : 22 Sep 21 02:58 pm using AcqMethod Hy_Pest_24cm_S10.M
Instrument : GCMS 8
Sample Name: E-018S
Misc Info : 2,4,6-Trichlorophenyl-4'-nitrophenyl ether; in iso
Vial Number: 61



AccuStandard®

Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5min, 10psi constant flow
Split 20:1, 1µl inj.; GC/MS; Injector 270°C

