

File :D:\MassHunter\GCMS\1\data\20230608\30- [JP].D
Operator : Organic QC Lab
Acquired : 08 Jun 2023 14:33 using AcqMethod Hy_Pest_15cm_S5.M
Instrument : GCMS 8
Sample Name: FRS-037S
Misc Info : FireMaster 680; in toluene
Vial Number: 35



Column: DB-5MS UI, 30m x 0.25mm ID x 0.25µm
Oven Program: 80°C, 17°C/min to 340°C (hold 8min)
GC Parameters: Pulsed split 40psi 0.5m in, 10psi constant flow
Split 2:1, 2µl inj.; GC/MS; Injector 270°C

